



# TOXIVEN BIOTECH PRIVATE LIMITED

Target: Sodium Channels

Format: Targeted Venom Discovery Array

Code: T-VDA<sup>Na+</sup>

## Product Description

The **sodium (Na<sup>+</sup>) channel Targeted Venom Discovery Array™ (T-VDA<sup>Na+</sup>)** is specifically designed to maximise discovery of new tools. Na<sup>+</sup> channels are important drug targets for a range of **neurological disorders**, specifically **pain**. Venoms from theraphosids (tarantulas), scorpions and snakes are rich sources of new Na<sup>+</sup> tools. These targeted arrays contain pure venom fractions from 12, 24, 48 or 96 species **optimised for identification of novel tools**. Each array contains characterised venoms active on sodium channels from the literature to act as positive controls. The control venoms for T-VDA<sup>Na+</sup> include *Thrixopelma puriens* (Peruvian velvet tarantula) where **Protox II**, a gating modifier of NaV1.7<sup>1</sup>, was discovered; *Androctonus australis* (Sahara scorpion) where several selective sodium channel tools have been discovered<sup>2</sup>; and *Crotalus durissus* (South American rattlesnake) venom which contains **crotamine**<sup>3</sup>, one of the very few snake-derived Na<sup>+</sup> channel toxins. The other venom fractions making up the library have been specially selected by our drug discovery scientists to maximise novel hit potential.

- Venoms are supplied lyophilised in Echo® qualified acoustic source plates (Labcyte Inc) and are useable on any SBS footprint liquid handling device or by hand.
- 384-well format has 200ng venom fraction per well, suggested dilution 20µl as hit fractions are typically active at 5µg/ml and below.
- 1536-well format also available.

1. Priest B.T., Blumenthal K.M., Smith J.J., Warren V.A., Smith M.M.(2007). ProTx-I and ProTx-II: gating modifiers of voltage-gated sodium channels. *Toxicon*, 49:194-201
2. Loret E.P., Martin-Eauclaire M.-F., Mansuelle P., Sampieri F., Granier C., Rochat H.(1991). An anti-insect toxin purified from the scorpion *Androctonus australis* hectoralso acts on the alpha-and beta-sites of the mammalian sodium channel: sequence and circular dichroism study. *Biochemistry* 30:633-640.
3. Mancin A.C., Soares A.M., Andriao-Escarso S.H., Faca V.M., Greene L.J., Zuccolotto S., Pela I.R., Giglio J.R.(1998). The analgesic activity of crotamine, a neurotoxin from *Crotalus durissus terrificus* (South American rattlesnake) venom: a biochemical and pharmacological study. *Toxicon*, 36:1927-1937

Data compiled from UniProt: Reorganizing the protein space at the Universal Protein Resource (UniProt), *Nucleic Acids Res.* 40: D71-D75 (2012).

### ToxiVen Biotech Private Limited

1-3B, First floor, MR Complex, Kovaipudur, Coimbatore, Tamil Nadu, India 641042

Tel: +91 422 436 7550 | Mob: +91 98940 70550

Email: [products@toxiven.com](mailto:products@toxiven.com) | Website: [www.toxiven.com](http://www.toxiven.com)